


Cgrp Rabbit Mab For Wb, Elisa - P06881

Item Number: CM0007307



Introduction

CamelBio's CGRP Rabbit mAb is a monoclonal antibody for WB and ELISA, targeting human Calcitonin gene-related peptide 1 (SWISS: P06881). CGRP is a vasodilatory hormone involved in cardiovascular and neuronal functions. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	CGRP Rabbit mAb
Remarks/Alias	CT; KC; PCT; CGRP; CALC1; CGRP1; CGRP-I; CGRP-alpha
Species	Human
Gene ID (human)	796
Gene ID	796
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 25-125 of human CGRP (NP_001029125.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P06881
Protein Weight	15kDa
Shipping	Ice bag