

Chdh Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q8Ne62

Item Number: CM0009076

Introduction

CHDH rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Targets mitochondrial choline dehydrogenase (65 kDa) for autophagy and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	CHDH Rabbit mAb
Remarks/Alias	CHDH
Species	Human
Gene ID (Human)	55349
Gene ID	55349
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 245-594 of human CHDH (NP_060867.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NE62
Protein Weight	65kDa
Shipping	Ice bag