

Chmp7 Rabbit Pab - Q8Wux9

Item Number: CM0017029

Introduction

Polyclonal CHMP7 antibody validated for WB, ELISA in human, mouse, rat. Targets ESCRT-III protein with nuclear envelope sealing function. Ideal for cell division and endosomal sorting studies.

[Learn More](#)

Parameter	Specification
Product Name	CHMP7 Rabbit pAb
Remarks/Alias	CHMP7
Species	Human
GeneID (Human)	91782
GeneID	91782
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WUX9
Protein Weight	51kDa
Shipping	Ice bag