

Chrna10 Rabbit Pab - Q9Gzz6

Item Number: CM0029362

Introduction

Rabbit polyclonal antibody against human CHRNA10, a neuronal nicotinic acetylcholine receptor subunit. Validated for Western blot and ELISA. Ideal for studying auditory synaptic transmission, keratinocyte adhesion, and calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	CHRNA10 Rabbit pAb
Remarks/Alias	CHRNA10
Species	Human
Gene ID (Human)	57053
Gene ID	57053
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-240 of human CHRNA10 (NP_065135.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZZ6
Protein Weight	50kDa
Shipping	Ice bag