


Cited2 Rabbit Polyclonal Antibody For Wb And Elisa - Q99967

Item Number: CM0020003



Introduction

CITED2 Rabbit Polyclonal Antibody validated for Western blotting and ELISA. Detects human CITED2 protein, cross-reacts with rat, predicted molecular weight 28kDa. Ideal for studying transcriptional regulation and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	CITED2 Rabbit pAb
Remarks/Alias	ASD8; MRG1; VSD2; MRG-1; P35SRJ; CITED2
Species	Human
GeneID (Human)	10370
GeneID	10370
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CITED2 (NP_001161860.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q99967
Protein Weight	28kDa
Shipping	Ice bag