

Citrulline-Histone H3 (R8 Citrullination) Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Db, Elisa

Item Number: CM0000078



Introduction

Rabbit monoclonal antibody targeting citrullinated Histone H3 at arginine 8 (R8). Validated for Western blot, dot blot, and ELISA. Cross-reacts with human and other species (wide range predicted). Ideal for epigenetic and citrullination studies.

[Learn More](#)

Parameter	Value
Product Name	Citrulline-Histone H3-R8 Rabbit mAb
Remarks/Alias	H3t, H3.4, H3/g, H3FT, H3C16, HIST3H3
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	A synthetic citrullinated peptide around R8 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag