

Claudin18.2 Rabbit Monoclonal Antibody For Ihc-P, Elisa - P56856

Item Number: CM0003728

Introduction

Specific to human CLDN18.2 isoform, this rabbit monoclonal antibody is validated for IHC-P and ELISA applications. Ideal for tight junction research, gastric and lung tissue studies, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Claudin18.2 Rabbit mAb
Remarks/Alias	SFTA5; SFTPJ; Claudin18.2
Species	Human
Gene ID	51208
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CLDN18.2 (P56856-2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P56856
Protein Weight	28kDa
Shipping	Ice bag