

Cldn18 Rabbit Pab - P56856

Item Number: CM0029971

Introduction

Reliable rabbit polyclonal antibody against human Claudin-18, validated for WB and ELISA. Shows cross-reactivity with mouse and rat, ideal for tight junction, lung, and stomach research.

[Learn More](#)

Parameter	Value
Product Name	CLDN18 Rabbit pAb
Remarks/Alias	SFTA5; SFTPJ; CLDN18
Species	Human
GeneID (Human)	51208
GeneID	51208
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 200-261 of human CLDN18 (NP_057453.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56856
Protein Weight	28kDa
Shipping	Ice bag