


Cldn2 Rabbit Polyclonal Antibody For Wb And Elisa - P57739

Item Number: CM0014065



Introduction

CLDN2 Rabbit Polyclonal Antibody validated for WB and ELISA in Human, Mouse, Rat. Targets Claudin-2, a tight junction protein involved in paracellular permeability. Ideal for epithelial barrier research.

[Learn More](#)

Parameter	Value
Product Name	CLDN2 Rabbit pAb
Remarks/Alias	OAZON; claudin-2; CLDN2
Species	Human
Gene ID (Human)	9075
Gene ID	9075
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57739
Protein Weight	25kDa
Shipping	Ice bag