

Cldn4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14493

Item Number: CM0014896



Introduction

Rabbit polyclonal antibody against human CLDN4, validated for Western blot, IHC paraffin, and ELISA. Recognizes Claudin-4, a tight junction protein involved in chloride permeability. Ideal for IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	CLDN4 Rabbit pAb
Remarks/Alias	CPER; CPE-R; CPETR; CPETR1; WBSCR8; hCPE-R; CLDN4
Species	Human
GeneID (Human)	1364
GeneID	1364
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 109-209 of human CLDN4 (NP_001296.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O14493
Protein Weight	22kDa
Shipping	Ice bag