


Cleaved Caspase-7 Rabbit Mab For Wb, Elisa - P55210

Item Number: CM0007174



Introduction

Rabbit monoclonal antibody against human cleaved Caspase-7, validated for Western blot and ELISA. Detects the active form of this executioner caspase in apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Cleaved Caspase-7 Rabbit mAb
Remarks/Alias	MCH3, CMH-1, LICE2, CASP-7, ICE-LAP3, Cleaved Caspase-7
Species	Human
GeneID (Human)	840
GeneID	840
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CASP7 (NP_001218.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P55210
Protein Weight	34kDa
Shipping	Ice bag