

Cleaved Gasdermin B (C Terminal) Rabbit Monoclonal Antibody For Wb And Elisa - Q8Tax9

Item Number: CM0007704



Introduction

Cleaved Gasdermin B (C Terminal) rabbit monoclonal antibody for Western blot and ELISA. Validated in human samples, targeting the C-terminal fragment after granzyme-mediated cleavage. Ideal for pyroptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Cleaved Gasdermin B (C Terminal) Rabbit mAb
Remarks/Alias	GSDML; PP4052; GSDMB-1; PRO2521; Cleaved Gasdermin B (C Terminal)
Species	Human
Gene ID (Human)	55876
Immunogen	Synthetic peptide corresponding to residues surrounding Ser245 of human Gasdermin B (NP_001159430.1)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TAX9
Protein Weight	18kDa/45kDa/46kDa/47kDa
Shipping	Ice bag