

## Cleaved Gasdermin B Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tax9-6

Item Number: CM0019435

### Introduction

Cleaved Gasdermin B Rabbit Polyclonal Antibody validated for WB and ELISA. Targets human Gasdermin B, a pore-forming protein mediating pyroptosis. Ideal for cell death and immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Cleaved Gasdermin B Rabbit pAb                                                                                                       |
| Remarks/Alias    | GSDML; PP4052; GSDMB-1; PRO2521; Cleaved Gasdermin B                                                                                 |
| Species          | Human                                                                                                                                |
| Gene ID          | 55876                                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human Cleaved Gasdermin B. (NP_001159431.1) |
| Source           | Rabbit                                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                                |
| Application      | WB, ELISA                                                                                                                            |
| Cross Reactivity | Human                                                                                                                                |
| SWISS            | Q8TAX9-6                                                                                                                             |
| Protein Weight   | 18kDa/45kDa/46kDa/47kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                                              |