


Cleaved Gasdermin D (C Terminal) Rabbit Mab - Q9D8T2

Item Number: CM0009577



Introduction

Cleaved Gasdermin D (C Terminal) Rabbit monoclonal antibody (Mouse) for WB, IP, ELISA. Targets the C-terminal region of gasdermin-D, a key pyroptosis mediator. Ideal for inflammasome and cell death research.

[Learn More](#)

Parameter	Value
Product Name	Cleaved Gasdermin D (C Terminal) Rabbit mAb
Remarks/Alias	DF5L; M2-4; Dfna5l; GsdmD-1; Gsdmdc1; 1810036L03Rik
Species	Mouse
Gene ID (Human)	79792
Gene ID	69146
Immunogen	Synthetic peptide corresponding to amino acids 277-376 of mouse Cleaved Gasdermin D (C Terminal) (NP_081236.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse
SWISS	Q9D8T2
Protein Weight	53 kDa
Shipping	Ice bag