


Cleaved Gasdermin D (Gly276) Rabbit Mab For Wb, Elisa - P57764

Item Number: CM0009338



Introduction

Cleaved Gasdermin D (Gly276) Rabbit mAb for WB and ELISA. Reacts with human Gasdermin D. Ideal for pyroptosis and inflammasome research. Monoclonal antibody from synthetic peptide immunogen.

[Learn More](#)

Parameter	Value
Product Name	Cleaved Gasdermin D Gly276 Rabbit mAb
Remarks/Alias	DF5L; DFNA5L; FKSG10; GSDMDC1
Species	Human
Gene ID (Human)	79792
Immunogen	Synthetic peptide A synthetic peptide corresponding to residues surrounding Asp276 of human Gasdermin D(NP_079012.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P57764
Protein Weight	53kDa
Shipping	Ice bag