


Cleaved IL-18 (Tyr37) Rabbit Mab For Wb - Q14116

Item Number: CM0007646



Introduction

Rabbit monoclonal antibody against cleaved human IL-18 at Tyr37, validated for Western blot. Detects a 22 kDa band. Ideal for pro-inflammatory cytokine research and pyroptosis studies.

[Learn More](#)

Parameter	Value
Product Name	Cleaved IL-18 (Tyr37) Rabbit mAb
Remarks/Alias	mature IL-18; Cleaved IL-18 (Tyr37)
Species	Human
GeneID (Human)	3606
GeneID	3606
Immunogen	Synthetic peptide corresponding to residues surrounding Tyr37 of human IL18 (NP_001553.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	Q14116
Protein Weight	22KDa
Shipping	Ice bag