

Cleaved Parp (Asp214) Rabbit Pab - Q921K2

Item Number: CM0018737

Introduction

Cleaved PARP (Asp214) Rabbit pAb detects the apoptosis-specific cleaved form of PARP in mouse samples. Validated for WB and ELISA, this polyclonal antibody targets amino acids 201-300.

[Learn More](#)

Product Name	Remarks/Alias	Species	Gene ID	Immunogen	Source	Category	Application	Cross Reactivity	SWISS	Protein Weight	Shipping
Cleaved PARP (Asp214) Rabbit pAb	PARP; PPOL; ARTD1; Adprp; Adprt1; msPARP; parp-1; sPARP-1; 5830444G22Rik; Cleaved PARP (Asp214)	Mouse	11545	A synthetic peptide corresponding to a sequence within amino acids 201-300 of Mouse PARP (Asp214) (NP_031441.2)	Rabbit	Polyclonal Antibodies	WB, ELISA	Mouse	Q921K2	113kDa	Ice bag