

Cleaved Parp1 P25 Rabbit Mab For Wb, Ip, Elisa - P09874

Item Number: CM0003838



Introduction

Cleaved PARP1 p25 Rabbit Monoclonal Antibody for Western Blot, Immunoprecipitation, and ELISA, with cross-reactivity to human and mouse. Useful for apoptosis and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Cleaved PARP1 p25 Rabbit mAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP; Cleaved PARP1 p25
Species	Human
Gene ID (Human)	142
Gene ID	142
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09874
Protein Weight	113 kDa
Shipping	Ice bag