


Cleaved Parp1 P25 Rabbit Polyclonal Antibody For Wb And Elisa - P09874

Item Number: CM0026600



Introduction

Rabbit polyclonal antibody against the cleaved p25 fragment of human PARP1, suitable for Western blot (WB) and ELISA. Specifically recognizes human PARP1, a key enzyme in DNA repair and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Cleaved PARP1 p25 Rabbit pAb
Remarks/Alias	PARP; PARS; PPOL; ADPRT; ARTD1; ADPRT1; PARP-1; ADPRT 1; pADPRT-1; Poly-PARP
Species	Human
Gene ID (Human)	142
Gene ID	142
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09874
Protein Weight	113kDa
Shipping	Ice bag