

Clec2B Rabbit Pab - Q92478

Item Number: CM0019977

Introduction

CLEC2B Rabbit Polyclonal Antibody (Q92478) validated for WB and ELISA. Detects human CLEC2B, cross-reacts with rat. Ideal for research into NK cell-mediated cytotoxicity and immune checkpoint regulation.

[Learn More](#)

Parameter	Value
Product Name	CLEC2B Rabbit pAb
Remarks/Alias	AICL; IFN1; CLECSF2; HP10085; CLEC2B
Species	Human
Gene ID (Human)	9976
Gene ID	9976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human CLEC2B (NP_005118.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q92478
Protein Weight	17kDa
Shipping	Ice bag