


## CLIC3 Rabbit Pab - O95833

Item Number: CM0017459



### Introduction

CLIC3 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, with reactivity against human, mouse, and rat CLIC3 (27kDa). It is suitable for studying chloride intracellular channel functions and tumor angiogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CLIC3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | CLIC3                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (human)   | 9022                                                                                                               |
| GeneID           | 9022                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 87-236 of human CLIC3 (NP_004660.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O95833                                                                                                             |
| Protein Weight   | 27kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |