

Cnap1 Rabbit Mab For Wb And Elisa - Q15021

Item Number: CM0012461

Introduction

High-quality CNAP1 (CAP-D2/hCAP-D2) rabbit monoclonal antibody validated for WB and ELISA. Detects human and mouse condensin complex subunit 1. Ideal for mitosis and chromosome condensation research. 157kDa, recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	CNAP1 Rabbit mAb
Remarks/Alias	CNAP1; CAP-D2; MCPH21; hCAP-D2
Species	Human
GeneID (Human)	9918
GeneID	9918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human CNAP1 (Q15021)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15021
Protein Weight	157kDa
Shipping	Ice bag