

Cntn4 Rabbit Mab For Wb, Ihc-P, Elisa - Q8Iwv2

Item Number: CM0011285

Introduction

CNTN4 Rabbit mAb is a rabbit monoclonal antibody for WB, IHC-P, and ELISA, detecting human, mouse, and rat CNTN4 (Contactin-4). It targets a neuronal cell adhesion molecule involved in neurite outgrowth and synapse formation.

[Learn More](#)

Parameter	Value
Product Name	CNTN4 Rabbit mAb
Remarks/Alias	AXCAM; BIG-2; CNTN4
Species	Human
Gene ID (Human)	152330
Gene ID	152330
Immunogen	Synthetic peptide corresponding to amino acids 800-900 of human CNTN4 (Q8I WV2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8I WV2
Protein Weight	113kDa
Shipping	Ice bag