

Cntn6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uq52

Item Number: CM0014679

Introduction

High-quality rabbit polyclonal antibody against human CNTN6 (Contactin-6, NB-3). Validated for Western blot and ELISA. Ideal for investigating nervous system development, oligodendrocyte generation, and Notch signaling. - Q9UQ52

[Learn More](#)

Parameter	Value
Product Name	CNTN6 Rabbit pAb
Remarks/Alias	NB3; CNTN6
Species	Human
GeneID (human)	27255
GeneID	27255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human CNTN6 (NP_055276.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UQ52
Protein Weight	114kDa
Shipping	Ice bag