

Coenzyme Nad Rabbit Monoclonal Antibody For Db And Elisa

Item Number: CM0000008

Introduction

This Coenzyme NAD Rabbit Monoclonal Antibody (A21047) is species-independent and suitable for DB and ELISA. It is designed for research use in NAD coenzyme detection.

[Learn More](#)

Parameter	Value
Product Name	Coenzyme NAD Rabbit mAb
Species	Species independent
GenelD	CAS:53-84-9
Immunogen	Chemical compounds NAD
Source	Rabbit
Category	Other Modified Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag