

Coenzyme Nad Rabbit Polyclonal Antibody For Db And Elisa - Nad

Item Number: CM0023501

Introduction

Rabbit polyclonal antibody targeting Coenzyme NAD, suitable for ELISA and dot blot (DB) applications. Species-independent reactivity enables broad NAD detection in diverse biological samples. Ideal for metabolism and redox research.

[Learn More](#)

Parameter	Value
Product Name	Coenzyme NAD Rabbit pAb
Immunogen	Chemical compounds NAD
Source	Rabbit
Category	Polyclonal Antibodies
Application	DB, ELISA
Cross Reactivity	Species independent
Shipping	Ice bag