


Col2A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02458

Item Number: CM0012828



Introduction

High-affinity COL2A1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for cartilage biology and skeletal development research.

[Learn More](#)

Parameter	Value
Product Name	COL2A1 Rabbit pAb
Remarks/Alias	AOM; ANFH; SEDC; STL1; COL11A3; COL2A1
Species	Human
Gene ID (Human)	1280
Gene ID	1280
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1220-1282 of human COL2A1 (NP_001835.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02458
Protein Weight	142 kDa
Shipping	Ice bag