

Col4A6 Rabbit Polyclonal Antibody For Wb, Elisa - Q14031

Item Number: CM0022536

Introduction

COL4A6 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human COL4A6, cross-reacts with rat. Suitable for extracellular matrix and basement membrane research.

[Learn More](#)

Parameter	Value
Product Name	COL4A6 Rabbit pAb
Remarks/Alias	DFNX6; DELXq22.3; CXDELq22.3; COL4A6
Species	Human
GeneID (Human)	1288
GeneID	1288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1575-1691 of human COL4A6 (NP_001838.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q14031
Protein Weight	164kDa
Shipping	Ice bag