

Col9A2 Rabbit Polyclonal Antibody For Wb, Elisa - Q14055

Item Number: CM0019012



Introduction

COL9A2 Rabbit Polyclonal Antibody (Q14055) raised against a synthetic peptide from human collagen alpha-2(IX) chain. Suitable for WB and ELISA with human and rat samples. Ideal for investigating cartilage, vitreous, and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	COL9A2 Rabbit pAb
Remarks/Alias	MED, EDM2, STL5, DJ39G22.4, COL9A2
Species	Human
Gene ID (Human)	1298
Gene ID	1298
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-689 of human COL9A2 (NP_001843.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q14055
Protein Weight	65kDa
Shipping	Ice bag