

Col9A3 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q14050

Item Number: CM0029381



Introduction

COL9A3 Rabbit pAb is a polyclonal antibody developed in rabbit against a recombinant human COL9A3 protein fragment. Validated for IHC-P and ELISA with mouse cross-reactivity, it supports collagen IX research in cartilage and ocular tissues.

[Learn More](#)

Parameter	Value
Product Name	COL9A3 Rabbit pAb
Remarks/Alias	IDD, MED, EDM3, STL6, DJ885L7.4.1, COL9A3
Species	Human
GeneID (human)	1299
GeneID	1299
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 445-684 of human COL9A3 (NP_001844.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q14050
Protein Weight	64kDa
Shipping	Ice bag