


Colgalt1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nbj5

Item Number: CM0019629



Introduction

COLGALT1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Targets procollagen galactosyltransferase 1 for collagen glycosylation and triple helix formation studies.

[Learn More](#)

Parameter	Value
Product Name	COLGALT1 Rabbit pAb
Remarks/Alias	B5VD3; GLT25D1; ColGalT 1; COLGALT1
Species	Human
Gene ID (Human)	79709
Gene ID	79709
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-230 of human COLGALT1 (NP_078932.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NBJ5
Protein Weight	72kDa
Shipping	Ice bag