


Collagen Iv Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q01955 / P29400 / P08572 / P02462

Item Number: CM0012451



Introduction

Rabbit monoclonal antibody targeting human Collagen IV (COL4A1,2,3,4) validated for WB, IF/ICC, ELISA. Ideal for basement membrane research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Collagen IV Rabbit mAb
Remarks/Alias	ATS2; ATS3; BFH2; ATS3A; ATS3B; BSVD; BSVD1; RATOR; PADMAL; COL4A1s
Species	Human
GeneID (Human)	1285/1282/1284/1286
GeneID	1285/1282/1284/1286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1555-1669 of human Collagen IV (NP_000082.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q01955/P29400/P08572/P02462
Protein Weight	162 kDa/161 kDa/168 kDa
Shipping	Ice bag