


Complement C3 Rabbit Mab For Wb, Elisa - P01027

Item Number: CM0005745



Introduction

Monoclonal rabbit antibody against mouse Complement C3 (P01027) for Western blot and ELISA. Cross-reacts with mouse and rat. Validated for research use in innate immunity and complement pathway studies.

[Learn More](#)

Parameter	Value
Product Name	Complement C3 Rabbit mAb
Remarks/Alias	ASP; Plp; HSE-MSF; Complement C3
Species	Mouse
GeneID (Human)	718
GeneID	12266
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01027
Protein Weight	186kDa
Shipping	Ice bag