

Complement C3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01024

Item Number: CM0028826



Introduction

Rabbit polyclonal antibody targeting human Complement C3, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for innate immunity and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Complement C3 Rabbit pAb
Remarks/Alias	ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p; Complement C3
Species	Human
Gene ID (Human)	718
Gene ID	718
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01024
Protein Weight	187kDa
Shipping	Ice bag