

Conjugated Aif1/Iba1 Rabbit Mab - O70200(Mouse),P55008(Human)

Item Number: CM0010568

Introduction

conjugated rabbit monoclonal antibody targeting AIF1/IBA1 (allograft inflammatory factor 1), validated for WB and IF-P applications with cross-reactivity to human, mouse, and rat. Suitable for immunoassays.

[Learn More](#)

Parameter	Value
Product Name	Conjugated AIF1/IBA1 Rabbit mAb
Remarks/Alias	G1, Iba1, AIF-1, D17H6S50E
Species	Mouse
Gene ID (Human)	199
Gene ID (Mouse)	11629
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse AIF1/IBA1 (O70200)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	O70200 (mouse), P55008 (Human)
Protein Weight	17kDa
Shipping	Ice bag