

Conjugated Aif1/Iba1 Rabbit Mab - O70200(Mouse),P55008(Human)

Item Number: CM0009613

Introduction

Conjugated AIF1/IBA1 Rabbit monoclonal antibody, validated for WB and IF-P. Detects AIF1/IBA1 in human, mouse, and rat. Ideal for neuroinflammation and microglia research, with consistent batch-to-batch performance for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Conjugated AIF1/IBA1 Rabbit mAb
Remarks/Alias	G1; Iba1; AIF-1; D17H6S50E
Species	Mouse
Gene ID (Human)	199
Gene ID (Mouse)	11629
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse AIF1/IBA1 (O70200)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	O70200(mouse), P55008(Human)
Protein Weight	17kDa
Shipping	Ice bag