

-Conjugated Anti-Cd200/Ox2 Monoclonal Antibody For Flow Cytometry - P41217

Item Number: CM0004887

Introduction

-conjugated rabbit monoclonal antibody targeting human CD200/OX2. Validated for flow cytometry, it specifically recognizes the cell surface glycoprotein involved in immune regulation. - P41217

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD200/OX2 mAb
Remarks/Alias	CD200; MOX1; MOX2; MRC; OX-2; CD200 molecule
Species	Human
Gene ID	4345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of human CD200/OX2 (NP_001352780.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P41217
Protein Weight	30kDa/31kDa/32kDa
Shipping	Ice bag