

-Conjugated Anti-Cd31/Pecam1 Monoclonal Antibody For Flow Cytometry - Q08481

Item Number: CM0007814

Introduction

-conjugated rabbit monoclonal antibody against mouse/rat CD31/PECAM1. Validated for flow cytometry, suitable for endothelial cell and leukocyte interaction studies. Targets PECAM1 adhesion molecule. Ideal for immunology and vascular biology research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
Gene ID (Human)	5175
Gene ID (Mouse)	18613
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	69kDa/77kDa/80kDa/81kDa
Shipping	Ice bag