

Conjugated Anti-Gfap Rabbit Monoclonal Antibody For If-P - P14136

Item Number: CM0008707



Introduction

Conjugated Anti-GFAP Rabbit Monoclonal Antibody validated for IF-P, reactive with human, mouse, and rat GFAP. Ideal for astrocyte marker studies in neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Conjugated GFAP Rabbit mAb
Remarks/Alias	ALXDRD
Species	Human
GeneID	2670
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GFAP (P14136))
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	P14136
Protein Weight	50kDa
Shipping	Ice bag