

Conjugated Tom20 Rabbit Mab For Wb, If/Icc, If-P - Q15388

Item Number: CM0005837



Introduction

-conjugated rabbit monoclonal antibody against human TOM20, validated for WB, IF/ICC, and IF-P. Recognizes human, mouse, and rat TOM20 for mitochondrial protein import and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	Conjugated TOM20 Rabbit mAb
Remarks/Alias	MAS20, MOM19, TOM20
Species	Human
GeneID (Human)	9804
GeneID	9804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-145 of human TOM20 (NP_055580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15388
Protein Weight	16kDa
Shipping	Ice bag