


Cop1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nhy2

Item Number: CM0023210



Introduction

COP1 Rabbit Polyclonal Antibody (pAb) validated for WB and ELISA. Suitable for investigating COP1's role in ubiquitination, p53 degradation, and cell survival pathways. Human cross reactivity.

[Learn More](#)

Parameter	Value
Product Name	COP1 Rabbit pAb
Remarks/Alias	FAP78; RFWD2; CFAP78; RNF200; COP1
Species	Human
Gene ID (Human)	64326
Gene ID	64326
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 582-731 of human COP1 (NP_071902.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NHY2
Protein Weight	80kDa
Shipping	Ice bag