

Copb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P53618

Item Number: CM0018628

Introduction

Rabbit polyclonal antibody raised against human COPB1 (Beta-COP). Validated for WB, IF/ICC, and ELISA; cross-reacts with human, mouse, rat. Ideal for studying coatomer-mediated Golgi/ER transport and lipid homeostasis.

[Learn More](#)

Parameter	Value
Product Name	COPB1 Rabbit pAb
Remarks/Alias	COPB; BARMACS; COPB1
Species	Human
GeneID (human)	1315
GeneID	1315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-740 of human COPB1 (NP_057535.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53618
Protein Weight	107kDa
Shipping	Ice bag