

Coq9 Rabbit Pab - O75208

Item Number: CM0017368



Introduction

Rabbit polyclonal antibody against human COQ9, a key player in ubiquinone biosynthesis. Validated for Western blot (WB) and ELISA applications, ensuring reliable detection of this 36 kDa mitochondrial target. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	COQ9 Rabbit pAb
Remarks/Alias	COQ10D5; C16orf49; COQ9
Species	Human
Gene ID (Human)	57017
Gene ID	57017
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human COQ9 (NP_064708.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75208
Protein Weight	36kDa
Shipping	Ice bag