


Cox Iv Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P13073

Item Number: CM0030854



Introduction

Validate human, mouse, and rat COX IV (COX4I1) with this rabbit polyclonal antibody, optimized for WB, IHC-P, IF/ICC, IP, and ELISA. Ideal for mitochondrial and oxidative phosphorylation studies.

[Learn More](#)

Parameter	Value
Product Name	COX IV Rabbit pAb
Remarks/Alias	COX4; COXIV; COX4-1; COXIV-1; MC4DN16; COX IV-1; COX IV
Species	Human
GeneID (Human)	1327
GeneID	1327
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-98 of human COX IV (NP_001852.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13073
Protein Weight	20kDa
Shipping	Ice bag