

Cox1 Rabbit Polyclonal Antibody - Cox1

Item Number: CM0025330



Introduction

Rabbit polyclonal antibody against COX1 (Cyclooxygenase-1), reactive with human, mouse, and rat. Suitable for WB, IF, IHC-P, and ELISA applications. Immunogen corresponds to mouse COX1 amino acids 415-514. Ideal for inflammation research.

[Learn More](#)

Parameter	Details
Product Name	COX1 Rabbit pAb
Alias	CoxI; COX1
Species	Mouse
Gene ID (Human)	4512
Gene ID (Mouse)	17708
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 415-514 of mouse COX1 (NP_904330.1)
Host Species	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	56 kDa
Shipping	Ice bag