

Cralbp Rabbit Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P12271

Item Number: CM0004769



Introduction

CRALBP rabbit monoclonal antibody for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, and rat CRALBP. Ideal for visual cycle and retinal research. 36 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	CRALBP Rabbit mAb
Remarks/Alias	CRALBP
Species	Human
GeneID (Human)	6017
GeneID	6017
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CRALBP (P12271).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12271
Protein Weight	36kDa
Shipping	Ice bag