

Crem Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q03060

Item Number: CM0027304

Introduction

High-quality rabbit polyclonal antibody against human CREM, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying spermatogenesis, circadian rhythm, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	CREM Rabbit pAb
Remarks/Alias	ICER; CREM-2; hCREM-2; CREM
Species	Human
GeneID (Human)	1390
GeneID	1390
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CREM (NP_853549.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03060
Protein Weight	37kDa
Shipping	Ice bag