

Crtap Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075718

Item Number: CM0026498



Introduction

High-quality rabbit polyclonal antibody targeting human CRTAP (cartilage-associated protein) for WB, IHC-P, and ELISA. Validated in human, mouse, and rat. Ideal for studying collagen hydroxylation and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	CRTAP Rabbit pAb
Remarks/Alias	OI7; CASP; P3H5; LEPREL3; CRTAP
Species	Human
Gene ID (Human)	10491
Gene ID	10491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human CRTAP (NP_006362.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75718
Protein Weight	47kDa
Shipping	Ice bag