

Csta Rabbit Pab - P01040

Item Number: CM0030195



Introduction

CSTA Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA, recognizing human and rat Cystatin-A (UniProt P01040). Suitable for protease inhibitor and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CSTA Rabbit pAb
Remarks/Alias	AREI; PSS4; STF1; STFA; CSTA
Species	Human
GeneID (human)	1475
GeneID	1475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human CSTA (NP_005204.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P01040
Protein Weight	11kDa
Shipping	Ice bag