

Ctbp2 Rabbit Pab For Wb, If/Icc, Elisa - P56545

Item Number: CM0021341



Introduction

CTBP2 Rabbit pAb is a polyclonal antibody raised against a C-terminal peptide of human CTBP2, suitable for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for transcription regulation and synapse research.

[Learn More](#)

Parameter	Value
Product Name	CTBP2 Rabbit pAb
Remarks/Alias	CTBP2
Species	Human
GeneID (Human)	1488
GeneID	1488
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human CTBP2 (NP_001320.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56545
Protein Weight	49kDa
Shipping	Ice bag